

3M

Cavilon™

Skin Care Solutions



3M™ Cavilon™ Advanced Skin Protectant

**Experience the power
of ultimate protection.**



The significance of healthcare-acquired skin damage.

Healthcare-acquired skin damage represents negative clinical outcomes resulting in potential complications such as infection, pain and suffering, and a poor patient experience. In addition, skin damage increases the work and cost of care.¹

Exposure over time to factors such as irritants, moisture, friction, shear and adhesives can lead to conditions of skin breakdown, including:

Moisture-Associated Skin Damage (MASD)



Incontinence-Associated Dermatitis (IAD)

Severe inflammation caused by liquid stool, mixed incontinence or urine, which can lead to destruction of the epidermis.



Peristomal Skin Damage

Problem stomas, poor stoma location and high-volume output, especially that of liquid stool, can contribute to skin injury that can rapidly progress to erosion.



Periwound Skin Damage

This type of skin damage is often associated with wounds that produce large quantities of drainage, such as venous ulcers or infected wounds.



Intertriginous Dermatitis (ITD)

Skin damage between skin surfaces due to the interaction of friction and moisture.

Pressure Ulcer/Injury



Localized damage to the skin and underlying soft tissue, usually over a bony prominence or related to a medical or other device. Moisture, friction and shear are accepted risk factors for pressure ulcer/injury development.

Medical Adhesive-Related Skin Injury (MARSI)




Skin damage, such as stripping or maceration, related to adhesive product use.

Limitations of traditional skin protectants.

Moisture barrier creams, ointments and pastes have long been the standard of care for skin protection, but often these products:

- Aren't effective for preventing and managing skin damage
- Don't last as long as you need them to
- Don't adhere to wet, weepy, damaged skin
- May interfere with healing
- Cause discomfort upon application, during wear and cleansing
- Don't stay in place on the skin
- Are difficult to clean and remove, and may cause additional skin damage
- Interfere with skin assessment
- Are not suitable around an ostomy or fistula



Newer technologies such as pure cyanoacrylates can have limitations as well. While able to attach to a wet surface, heat is given off during drying and the surface can fracture easily, compromising protection and requiring reapplication.

The power of 3M™ Cavilon™ Advanced Skin Protectant.

Featuring 3M's revolutionary polymer-cyanoacrylate technology, 3M™ Cavilon™ Advanced Skin Protectant is designed to help manage moderate to severe skin damage and protect at-risk skin. The ultra-thin yet highly durable barrier is able to attach to wet, weepy surfaces and create a protective environment that repels irritants and supports healing.

Attaches to wet, weepy, damaged skin¹

Provides an effective barrier, which has been shown to reduce the pain of managing IAD²

Single-use applicator reduces the potential for cross-contamination

Lasts up to 7 days³

3M™ Cavilon™ Advanced Skin Protectant is available in two sizes:

- 2.7ml applicator coverage area of approx. 100 sq. inches
- 0.7ml applicator coverage area of approx. 25 sq. inches
 - Ideal for application around ostomies, tubes and wounds

The results are clear

A clear improvement on traditional treatment options, 3M™ Cavilon™ Advanced Skin Protectant uses a unique polymer-cyanoacrylate system to deliver ultimate protection and prevention – even in the most challenging circumstances.



Patient with skin damage on Day 2



Patient with improvement noted on Day 4



*Product is not an analgesic.

See the science of skin protection at work.

3M™ Cavilon™ Advanced Skin Protectant represents a revolutionary technology for management of skin damage and protection of at-risk skin. Its formulation is unlike any other skin protectant or moisture barrier, but what makes it so different?



Unique, elastomeric polymer

The polymer forms a coating that can elongate and conform, avoiding the cracking seen with other moisture barriers. This assures greater barrier integrity, durability and protection against challenging irritants such as liquid stool and gastric fluid.



Revolutionary polymer-cyanoacrylate system

The cyanoacrylate enables attachment of the skin protectant to damaged skin that is wet and weeping. Once on the skin, the protective coating creates an environment that repels irritants and supports healing and comfort.

Remember: if a moisture barrier product cannot reliably attach to the underlying skin, it's not capable of reliable protection.




Non-stinging solvent

The polymer-cyanoacrylate system is delivered onto the skin by a non-stinging solvent.

3M™ Cavilon™ Advanced Skin Protectant is also breathable, allowing for moisture-vapor transmission that helps keep skin comfortable. Plus, it doesn't require removal, and the surface is easily cleansed – making wear easier for patients and care easier for staff.

How 3M™ Cavilon™ Advanced Skin Protectant works





The many uses of 3M™ Cavilon™ Advanced Skin Protectant.

Moisture-Associated Skin Damage (MASD)

Pressure Ulcer/Injury

Medical Adhesive-Related Skin Injury (MARS)



Incontinence-Associated Dermatitis (IAD)

Peristomal Skin Damage

Periwound Skin Damage

Intertriginous Dermatitis (ITD)

Manage damaged or broken skin

Manage IAD	Manage peristomal/perifistular skin damage	Manage periwound skin damage (e.g. maceration)	Manage superficial skin damage from moisture and friction	Manage superficial skin injury in difficult-to-dress locations	Manage superficial skin damage (e.g. stripping, skin tears) from adhesive use
------------	--	--	---	--	---

Protect at-risk skin

Protect intact skin especially in the presence of diarrhea or mixed incontinence	Protect skin around problem fecal or urinary stomas, fistulas or tracheostomies	Protect skin around at-risk wounds (e.g. heavily draining wounds such as diabetic foot ulcers, venous leg ulcers or infected wounds)		Protect intact skin from moisture, friction or shear	
--	---	--	--	--	--

Coverage area comparison:



3M™ Cavilon™ Advanced Skin Protectant 5050
2.7ml = 100 sq. in.
 applicator approx. coverage



3M™ Cavilon™ Advanced Skin Protectant 5051
0.7ml = 25 sq. in.
 applicator approx. coverage

Ultimate protection that can effectively transform your practice.



For your patients and residents, 3M™ Cavilon™ Advanced Skin Protectant provides an effective barrier, which can help:

- ✓ Prevent skin injury and maintain skin integrity
- ✓ Protect against caustic, corrosive body fluids including liquid stool and gastric fluid
- ✓ Create an environment that supports healing
- ✓ Improve the overall healthcare experience



For you and your facility, 3M™ Cavilon™ Advanced Skin Protectant can help:

- ✓ Eliminate the need for frequent reapplication
- ✓ Reduce the potential for the cross-contamination that's possible with traditional multi-use products
- ✓ Offer superior value in dollars spent for skin management
- ✓ Lead to improved patient and resident outcomes



100% of patients

who reported pain associated with IAD on Day 1 (n=9) experienced a reduction in pain resulting from the protective barrier provided by 3M™ Cavilon™ Advanced Skin Protectant.*

*Product is not an analgesic.

With 3M™ Cavilon™ Advanced Skin Protectant, your facility could potentially save (per IAD patient per week)*:



\$533



16.5 hours of nursing time

*Based on six episodes per day⁴, two nurses at 3.9 hours/week at \$33.23/hour⁵ each; application/cleansing time of 45 sec/2 min⁶; three applications per week.

Experience the power of
3M™ Cavilon™ Advanced Skin Protectant
at [3M.com/CavilonAdvanced](https://www.3M.com/CavilonAdvanced)

Ordering Information

3M™ Cavilon™ Advanced Skin Protectant

Cat. No	Size	Approx. Coverage	Ideal for:	Items/Box
5050	2.7ml applicator	100 sq. in.	<ul style="list-style-type: none">Protecting skin around fistulaIAD prevention and management	20
5051	0.7ml applicator	25 sq. in.	<ul style="list-style-type: none">Peristomal, peri tube and periwound protection	20

Purchasing 3M™ Cavilon™ Products

To order Cavilon products, visit 3M.com/Cavilon to find a distributor, contact your 3M Critical & Chronic Care Representative or call the 3M Helpline at 1-800-228-3957. Outside of the United States, contact the local 3M subsidiary.



Discover all the ways 3M™ Cavilon™ Skin Care Solutions can help you transform patient skin integrity at 3M.com/Cavilon

¹Brennan, Mary R.; Milne, Catherine T.; Agrell-Kann, Marie; Ekholm, Bruce P. Clinical Evaluation of a Skin Protectant for the Management of Incontinence Associated Dermatitis: An Open-Label, Nonrandomized, Prospective Study. *J of Wound, Ostomy & Continence Nursing*. 2017. 44(2):172-180.

²Been RA et al. *In vivo* methods to evaluate a new skin protectant for loss of skin integrity. *Wound Repair and Regen*. 2016; 24: 851-859.

³3M data on file. EM-05-01 3924.

⁴Bliss DZ, Zehrer C, Savik K, Smith G, Hedblom E. An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence. *J WOCN* 2007;34(2):143-52.

⁵Bureau of Labor Statistics. Occupational Employment Statistics for 31-1014 Nursing Assistant. (May 2015). Accessed on 5/18/2016 <http://www.bls.gov/oes/current/oes311014.htm>

⁶Heidegger CP; Graf S; Perneger T; Genton L; Oshima T; Pichard C. The burden of diarrhea in the intensive care unit (ICU-BD). A survey and observational study of the caregivers' opinions and workload. *Int J Nurs Stud*. 2016 Jul;59:163-8.



Medical Solutions Division

3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

3M and Cavilon are registered trademarks of 3M.
Please recycle. Printed in USA © 3M 2019.
All rights reserved. 70-2011-7076-1