

One viable tissue

- + Two preservation methods
- + Three natural components

Placental membranes for wound management composed of native viable cells, growth factors, and extracellular matrix

**Smith+Nephew**

**GrafixPL** <sup>◇</sup>

Lyopreserved Placental Membrane

**Grafix** <sup>◇</sup>

Cryopreserved Placental Membrane

# + Did you know?

Wounds can be devastating to your patients' lives and hospital outcomes. Major disability, poor quality of life, and reduced productivity contribute to the growing financial burden.<sup>1-3</sup> When healing stalls, so do the lives of your patients.

## Diabetic foot ulcers

1.5 

million U.S. patients affected each year<sup>4</sup>

85% 

of lower limb amputations in patients with diabetes are preceded by ulcerations<sup>5-7</sup>

40% 

recurrence rate within 1 year<sup>8</sup>

## Venous leg ulcers

1 

million U.S. patients affected each year<sup>4</sup>

66% 

fail to heal with standard of care (SOC) in 12 weeks<sup>9-11</sup>

70% 


recurrence rate within 3 months<sup>12</sup>

Wound duration is an independent **risk factor** for infection, hospitalization and amputation, contributing to significant patient morbidity and cost of care




# + What if there was a better solution for your patient?


An evidence-based approach for wound preparation to advance chronic wounds toward healing<sup>13-15</sup>




Tissue is non-viable or deficient  
Debridement<sup>13-16</sup>



Infection, inflammation, and biofilm  
Bioburden management<sup>13-15</sup>




Moisture imbalance  
Exudate management<sup>13-15</sup>



**Epithelial edge advancement**  
Promote epithelialization<sup>13-15</sup>

## Placental membranes may promote re-epithelialization of chronic wounds<sup>17</sup>

- Advanced cryopreservation and lyopreservation methods preserve the native properties and components of fresh placental membranes
- Other tissue preserving methods may alter or destroy the components and/or properties of fresh placental membranes<sup>18-21</sup>



\*List not inclusive of all growth factors.



†Present in amnion only.

# GrafixPL<sup>◇</sup> + Grafix<sup>◇</sup>

## About

- Lyopreserved or cryopreserved placental membranes derived from the amnion or chorion placental membranes
  - GrafixPL is lyopreserved and stored at room temperature
  - Grafix is cryopreserved and stored at -75°C to -85°C
- Trophoblast layer and maternal components removed to prevent an immune response
- Can be used as a wound cover or surgical wrap

## Cross-section of the placental membranes




Hydrated GrafixPL PRIME and thawed Grafix PRIME are equivalent<sup>22</sup>

# Engineered by nature + optimally preserved + available on demand

Many placental membrane products are available, but most have only growth factors and an extracellular matrix. All three native components of placental membranes, including viable cells, are preserved in GrafixPL<sup>o</sup> and Grafix<sup>o</sup>.

Endogenous cells remain viable<sup>22</sup>

3D matrix remains intact<sup>22</sup>



**LIVE** and **DEAD** cell staining

The presence of viable cells in GrafixPL and Grafix was independently confirmed by researchers at Rutgers University, Montana State University and the University of Texas Southwestern.<sup>23-24</sup>

Native growth factors are retained<sup>22</sup>

	Fresh amnion	GrafixPL PRIME	Grafix PRIME
IL-10	■	■	■
IL-1RA	■	■	■
PDGF-BB	■	■	■
bFGF	■	■	■
SDF-1α	■	■	■
Angiopoietin-1	■	■	■

# Quality evidence + Results driven

Grafix<sup>◇</sup> delivered statistically significant improvements over standard wound care for closing diabetic foot ulcers<sup>25</sup>

## Study overview<sup>25</sup>

- Prospective, 20-center, randomized, single-blinded, controlled trial\*
- Third party blinded image verification

## Results<sup>25</sup>



Significantly higher complete closure rate

**(62% vs. 21%)**

- Faster median time to complete closure (42 days vs. 70 days)
- Decreased number of treatments (6 vs. 12)\*
- Fewer wound complications (44% vs. 66%)
- Fewer wound-related infections (18% vs. 36%)
- **65% complete closure** in the open-label crossover phase for patients who previously failed with standard treatment in the control group<sup>26</sup>
- Fewer infection-related hospitalizations (6% vs. 15%)

## Consideration<sup>25</sup>

- At the pre-specified interim analysis, study enrollment was terminated at the recommendation of the blinded review committee due to the superiority of Grafix versus the control



\*Trial had three phases: 1) a blinded phase of treatment with weekly Grafix (n=50) vs treatment with standard of care alone (n=47) with a primary endpoint of 100% re-epithelialization by week 12; 2) a follow-up phase with visits every 4 weeks for an additional 12 weeks; and 3) an open-label phase where patients in the control arm failing to close after 12 weeks were given the opportunity to receive Grafix for 12 weeks.

# Faster wound closure

## + Lower costs<sup>25,27\*</sup>



### \$14,813

lower cost of care for Grafix<sup>®</sup>-treated patients (n=50) compared to control (n=47, standard of care alone)

### \$13,828

lower cost of care for closed wounds (n=41) compared to non-closed wounds (n=56)

The lower costs for Grafix-treated patients were driven by faster wound closure, fewer wound complications, and fewer hospitalizations

## Low bias

## + More confidence<sup>28</sup>

### Assessment overview

The National Institute for Health and Care Excellence (NICE) provides national guidance and advice to improve health and social care. NICE evaluated the quality of data for dermal substitute studies that included cure rates at 12 weeks vs standard care in diabetic foot ulcers.

### Results

Grafix had the **highest overall effect (z-score)** and was the only study rated as **high quality** with no serious risk of bias

Product	z-score <sup>29</sup>	p-value <sup>29</sup>	NICE quality of evidence rating <sup>30</sup>
Grafix <sup>25</sup>	z=3.55	p=0.0004	High quality with low risk of bias
Dermagraft <sup>31-33</sup>	z=3.13	p=0.002	Moderate quality with serious risk of bias
EpiFix <sup>34</sup>	z=2.13	p=0.03	Low quality with very serious risk of bias
Apligraf <sup>35</sup>	z=2.10	p=0.04	Very low quality with serious risk of bias

\*Comparing estimated healthcare costs associated with Grafix versus standard of care alone in the RCT reported by Lavery et al. The cost of care was calculated based on treatments, medications, clinical procedures, and wound complications.



# Results in a large, retrospective, WoundExpert analysis of Grafix<sup>◇</sup> in DFUs mirror RCT closure rates<sup>36</sup>

## Study overview<sup>36</sup>

- Retrospective, 58-center analysis of Grafix in the management of DFUs with Net Health's WoundExpert electronic health record (EHR) database
- Population: All patients who received Grafix in the management of DFUs over a 4-year period were evaluated\* (360 patients, 441 wounds)

## Patient demographics and wound characteristics

Multiple wounds <sup>†</sup> in 4-year period	Approx. 90%
>6 wounds <sup>†</sup> in 4-year period	Approx. 50%
Wound size (mean)	5.1 cm <sup>2</sup>
Wound duration prior to treatment (mean)	102.4 days
Complex wounds with exposed bone, tendon, or joint capsule	14.7%

## Results<sup>36,‡</sup>



Complete wound closure (end of treatment)

**59.4%**



Time to wound closure (mean)

**57.7 days**



Number of grafts to close (mean)

**5.2**



Safety outcomes:

Amputations

**3%**

Related infections

**2%**

## Comparison between EHR real world study and randomized, controlled trial

Study type	Retrospective, multicenter <sup>36</sup>	Prospective, multicenter RCT <sup>25</sup>
Wounds	350	97 (50 Grafix, 47 Control)
Complex wounds	Allowed	Excluded
Complete wound closure at end of treatment	<b>59.4%</b>	<b>62.0%</b>

\*Exclusion criteria: Wounds missing baseline/follow-up measurements or receiving other skin substitute treatment concurrent with Grafix.

†DFUs and wounds of other etiologies.

‡350 wounds out of the 441 wounds evaluated were analyzed for closure and closure-related outcomes. Wounds  $\leq 0.25$  cm<sup>2</sup> were not included in the closure analysis. 0.25 cm<sup>2</sup> was chosen as a minimum size for closure analysis since in previous publications of WoundExpert database studies, wounds  $\leq 0.25$  cm<sup>2</sup> were considered closed.<sup>37-39</sup>



# Overwhelming success in the management of chronic complex DFUs<sup>40</sup>

## Study overview<sup>40</sup>

- Prospective, multicenter, open-label, single-arm trial of Grafix<sup>®</sup> in the management of DFUs with exposed bone and/or tendon
- 31 patients enrolled, 27 patients completed\*

## Wound characteristics

Wound duration prior to study (mean)	7.5 months
Wound size (mean)	14.6 cm <sup>2</sup>
Prior advanced wound therapy	67.7%

## Results<sup>40</sup>



Complete granulation at 16 weeks:

# 96.3%

Mean 6.8 graft applications in 6.8 weeks to achieve 100% granulation


Complete wound closure at 16 weeks:

# 59.3%

Mean 9.0 graft applications in 9.1 weeks to achieve complete wound closure

### 100% Granulation

Case example 1




Week 0 wound area: 70 cm<sup>2</sup>

Week 16 wound area: 0.8 cm<sup>2</sup>

### 100% Re-epithelialization


Case example 3



Week 0 wound area: 17.4 cm<sup>2</sup>

Week 12 wound closure


Case example 2



Week 0 wound area: 47.2 cm<sup>2</sup>

Week 16 wound area: 1.2 cm<sup>2</sup>

Case example 4



Week 0 wound area: 10.8 cm<sup>2</sup>

Week 10 wound closure

\*Two patients withdrew for non-compliance and two for surgical intervention.

This series of studies is provided for informational and educational purposes only. These cases may not represent typical outcomes. Every procedure and each patient undergoing wound care treatment represents unique sets of circumstances and, therefore, results may vary. Smith & Nephew does not provide medical advice. The information presented is not, and is not intended to serve as, medical advice. It is the responsibility of the treating physician to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients.

# Right product + Right time

## Grafix<sup>◇</sup> helped improve closure of chronic refractory venous leg ulcers<sup>41</sup>

### Study overview<sup>41</sup>

- Prospective longitudinal crossover study of Grafix in the management of refractory chronic VLU
- Single-center, open-label, single-arm where each patient served as their own control



81 Patients entered 12 week standard therapy phase\*

21 patients with 30 VLUs failed to heal with 12 weeks of standard therapy

They crossed over to the 12 week Grafix treatment phase<sup>†</sup>

### Results<sup>41</sup>

VLUs that failed to close with 12 weeks of standard therapy made **significantly greater progress** toward closure when **Grafix** was added to the treatment

Complete wound closure of chronic refractory VLUs:

**53%**

Mean 7.2 graft applications in

**10.9 weeks**

12 week follow-up phase

**no recurrence**

#### Outcomes of standard therapy phase vs Grafix treatment phase in 21 crossover patients

	Standard therapy phase	Grafix treatment phase	p-value
Baseline wound size (mean)	17.1 cm <sup>2</sup>	12.2 cm <sup>2</sup>	
Complete wound closure	0%	53%	<0.001
Wound area reduction (mean)	29%	79%	<0.001

\*Standard therapy phase: Patients were treated for 12 weeks with SOC (included multi-layer compression).

<sup>†</sup>Inclusion criteria for the Grafix treatment phase included: 1) failed to heal in standard therapy phase; 2) venous insufficiency confirmed by duplex ultrasound; 3) no infection, ischemia, or immunosuppression; and 4) radiofrequency ablation of the great saphenous vein for patients with evidence of superficial venous insufficiency. Radiofrequency ablation of the ipsilateral great saphenous vein was performed in 14 of the 21 patients at 4 weeks (mean) prior to entering the study.

# Chronic wound closure outcomes with GrafixPL<sup>◇</sup> suggest clinical equivalency with Grafix<sup>◇42</sup>

## Study overview<sup>42</sup>

- Retrospective, open-label, 5-center study on GrafixPL PRIME in the management of chronic wounds\*
- Wounds included DFUs (n=41), VLU's (n=19), surgical wounds (n=10), & other wounds (n=28)<sup>†</sup>

## Patient demographics and wound characteristics

Patients	78
Wounds	98
Wound size (mean)	13.3 cm <sup>2</sup>
Wound duration prior to treatment (mean)	8.7 months



## Results<sup>42</sup>

Complete wound closure<sup>42</sup>

**59.2%**

Median 6 graft applications in  
**63 days**  
to closure

**GrafixPL closure rates similar to those previously reported for Grafix**



\*Defined as wounds with no progression toward closure with 4 weeks of SOC or wounds in patients with significant comorbidities that put them at high risk for nonclosure.

<sup>†</sup>Other wounds include pressure ulcers, arterial wounds, chronic wounds, open hematomas, gangrenous wounds, radiation necrosis wounds, lymphedema wounds, ischemic wounds, & necrotizing fasciitis.


# GrafixPL<sup>◊</sup> and Grafix<sup>◊</sup> are placental membranes composed of native living cells, growth factors, and an intact extracellular matrix

- Designed for application directly to wounds and compromised surgical sites
- Flexible, conforming covers that may be applied over bone, tendon, and other structures
- GrafixPL and Grafix are available in multiple sizes, helping you reduce waste and cost

Placental tissue source	Product description	Part #
Lyopreserved amniotic membrane	GrafixPL PRIME 16 mm Disc (2 cm <sup>2</sup> )	PS13016
	GrafixPL PRIME 1.5 cm x 2 cm (3 cm <sup>2</sup> )	PS13015
	GrafixPL PRIME 2 cm x 3 cm (6 cm <sup>2</sup> )	PS13023
	GrafixPL PRIME 3 cm x 3 cm (9 cm <sup>2</sup> )	PS13033
	GrafixPL PRIME 3 cm x 4 cm (12 cm <sup>2</sup> )	PS13034
	GrafixPL PRIME 5 cm x 5 cm (25 cm <sup>2</sup> )	PS13055

Placental tissue source	Product description	Part #
Cryopreserved amniotic membrane	Grafix PRIME 16 mm Disc (2 cm <sup>2</sup> )	PS60013
	Grafix PRIME 1.5 cm x 2 cm (3 cm <sup>2</sup> )	PS11015
	Grafix PRIME 2 cm x 3 cm (6 cm <sup>2</sup> )	PS11023
	Grafix PRIME 3 cm x 4 cm (12 cm <sup>2</sup> )	PS11034
	Grafix PRIME 5 cm x 5 cm (25 cm <sup>2</sup> )	PS11055
Cryopreserved chorionic membrane	Grafix CORE 3 cm x 4 cm (12 cm <sup>2</sup> )	PS12034
	Grafix CORE 5 cm x 5 cm (25 cm <sup>2</sup> )	PS12055

## Actual sizes shown



References: 1. Sen CK, Gordillo GM, Roy S, et al. Human skin wounds: a major and snowballing threat to public health and the economy. *Wound Repair Regen.* 2009; 17(6): 763-771. 2. Fife AC, Carter MJ, Walker D, et al. Wound care outcomes and associated cost among patients treated in the US outpatient wound centers: Data from the US Wound Registry. *Wounds.* 2012; 24(1): 10-17. 3. Greer N, Foman NA, MacDonald R, et al. Advanced wound care therapies for nonhealing diabetic, venous, and arterial ulcers: a systematic review. *Ann Intern Med.* 2013; 159(8): 532-42. 4. SMARTTrak™. <https://app.smarttrak.com/markets/qs/2f491c456352f41ac0a80237b1942c19>. Retrieved 07 Nov 2019. 5. Frykberg RG, Zgonis T, Armstrong DG, et al. Diabetic foot disorders. A clinical practice guideline (2006 revision). *J Foot Ankle Surg.* 2006; 45(5): S1-S66. 6. Singh N, Armstrong DG, Lipsky BA. Preventing foot ulcers in patients with diabetes. *JAMA.* 2005; 293(2): 217-228. 7. Apelqvist J, Larsson J. What is the most effective way to reduce incidence of amputation in the diabetic foot? *Diabetes Metab Res Rev.* 2000; 16(S1): S75-S83. 8. Armstrong DG, Boulton AJM, Bus SA. Diabetic Foot Ulcers and Their Recurrence. *N Engl J Med.* 2017; 376(24): 2367-2375. 9. Mostow EN, Haraway GD, Dalsing M, et al. Effectiveness of an extracellular matrix graft (OASIS Wound Matrix) in the treatment of chronic leg ulcers: a randomized clinical trial. *J Vasc Surg.* 2005; 41(5): 837-43. 10. Harding K, Sumner M, Cardinal M. A prospective, multicentre, randomised controlled study of human fibroblast-derived dermal substitute (Dermagraft) in patients with venous leg ulcers. *Int Wound J.* 2013; 10(2): 132-7. 11. Bianchi C, Cazzell S, Vayser D, et al. A multicentre randomised controlled trial evaluating the efficacy of dehydrated human amnion/chorion membrane (EpiFix®) allograft for the treatment of venous leg ulcers. *Int Wound J.* 2018; 15(1): 114-122. 12. Franks PJ, Barker J, Collier M, et al. Management of Patients With Venous Leg Ulcers: Challenges and Current Best Practice. *J Wound Care.* 2016; 25 Suppl 6: S1-S67. 13. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. *Wound Repair Regen.* 2003; 11(Suppl 1): S1-S28. 14. Leaper DJ, Schultz G, Carville K, et al. Extending the TIME concept: what have we learned in the past 10 years? *Int Wound J.* 2012; 9(Suppl 2): 1-19. 15. Edmonds M, Foster AVM, Vowden P. Wound bed preparation for diabetic foot ulcers. In: *Position Document: Wound Bed Preparation in Practice.* London, England: MEP Ltd; 2004:6-11. 16. Kirsner RS. The standard of care for evaluation and treatment of diabetic foot ulcers. Available at: <https://barrystorage.blob.core.windows.net/assets/docs/continuing-medical-education/diabetic.pdf>. Accessed October 20, 2017. 17. Mamede AC, Carvalho MJ, Abrantes AM, et al. Amniotic membrane: from structure and functions to clinical applications. *Cell Tissue Res.* 2012; 349(2): 447-458. 18. Rodriguez-Ares MT, López-Valladares MJ, Touriño R, et al. Effects of lyophilization on human amniotic membrane. *Acta Ophthalmol.* 2009; 87(4): 396-403. 19. Lim LS, et al. Biological and ultrastructural properties of acelagraft, a freeze-dried-irradiated human amniotic membrane. *Arch Ophthalmol.* 2010; 128(10): 1303-1310. 20. Niknejad H, Deihim T, Solati-Hashjin M, et al. The effects of preservation procedures on amniotic membrane's ability to serve as a substrate for cultivation of endothelial cells. *Cryobiology.* 2011; 63(3): 145-151. 21. Thomasen H, Paulkin M, Steuhl KP, et al. Comparison of cryopreserved and air-dried human amniotic membrane for ophthalmologic applications. *Graefes Arch Clin Exp Ophthalmol.* 2009; 247(12): 1691-1700. 22. Dhall S, Sathyamoorthy M, Kuang JQ, et al. Properties of viable lyopreserved amnion are equivalent to viable cryopreserved amnion with the convenience of ambient storage. *PLoS One.* 2018; 13(10): e0204060. 23. Mao Y, Hoffman T, Dhall S, et al. Endogenous viable cells in lyopreserved amnion retain differentiation potential and anti-fibrotic activity in vitro. *Acta Biomater.* 2019; 94: 330-339. 24. Data on file at Osiris Therapeutics, Inc. 25. Lavery LA, Fulmer J, Shebetka KA, et al. The efficacy and safety of Grafix® for the treatment of chronic diabetic foot ulcers: results of a multi-centre, controlled, randomised, blinded, clinical trial. *Int Wound J.* 2014; 11(5): 554-560. 26. Lavery LA, Fulmer J, Shebetka KA, Regulski M, Vayser D, Fried D, Kashefsky H, Owings TM, Nadarajah J, Hesp Z. Open-label Extension Phase of a Chronic Diabetic Foot Ulcer Multicenter, Controlled, Randomized Clinical Trial Using Cryopreserved Placental Membrane. *Wounds Epub* 2018 June 29. 27. Nuccio EJ, Lavery LA, Min SJ. Innovative Treatment of Chronic Diabetic Foot Ulcer in a Controlled Randomized Clinical Trial Produces Fewer Adverse Events, Faster Wound Closure, and Lower Costs. *Diabetes Manag.* 2016; 6(4): 76-81. 28. Diabetic foot problems: prevention and management. National Institute for Health and Care Excellence (NICE) guideline. Published: 26 August 2015. <https://www.nice.org.uk/guidance/ng19>. 29. Diabetic foot problems: prevention and management. National Institute for Health and Care Excellence (NICE) guideline. Published: January 2016. <https://www.nice.org.uk/guidance/ng19/evidence/appendix-h-data-analysis-435358773>, page 11. 30. Diabetic foot problems: prevention and management. National Institute for Health and Care Excellence (NICE) guideline. Published: January 2016. <https://www.nice.org.uk/guidance/ng19/evidence/appendix-i-grade-profiles-15672915540>. 31. Gentzkow GD, Iwasaki SD, Hershon KS, et al. Use of dermagraft, a cultured human dermis, to treat diabetic foot ulcers. *Diabetes Care.* 1996; 19(4): 350-4. 32. Hanft JR, Surprenant MS. Healing of chronic foot ulcers in diabetic patients treated with a human fibroblast-derived dermis. *J Foot Ankle Surg.* 2002; 41(5): 291-9. 33. Marston WA, Hanft J, Norwood P, et al. The efficacy and safety of Dermagraft in improving the healing of chronic diabetic foot ulcers: results of a prospective randomized trial. *Diabetes Care.* 2003; 26(6): 1701-5. 34. Zelen CM, Serena TE, Denoziero G, et al. A prospective randomised comparative parallel study of amniotic membrane wound graft in the management of diabetic foot ulcers. *Int Wound J.* 2013; 10(5): 502-7. 35. Edmonds M, European and Australian Apligraf Diabetic Foot Ulcer Study Group. Apligraf in the treatment of neuropathic diabetic foot ulcers. *Int J Low Extrem Wounds.* 2009; 8(1): 11-8. 36. Raspovic KM, Wukich DK, Naiman DQ, Lavery LA, Kirsner RS, Attinger CE, Danilkovitch A. Effectiveness of viable cryopreserved placental membranes for management of diabetic foot ulcers in a real world setting. *Wound Repair Regen.* 2018 Apr 23. doi: 10.1111/wrr.12635. 37. Kirsner RS, Sabolinski ML, Parsons NB, Skornicki M, Marston WA. Comparative effectiveness of a bioengineered living cellular construct vs. a dehydrated human amniotic membrane allograft for the treatment of diabetic foot ulcers in a real world setting. *Wound Repair Regen.* 2015; 23(5): 737-44. 38. Marston WA, Sabolinski ML, Parsons NB, Kirsner RS. Comparative effectiveness of a bilayered living cellular construct and a porcine collagen wound dressing in the treatment of venous leg ulcers. *Wound Repair Regen.* 2014; 22(3): 334-40. 39. Kraus I, Sabolinski ML, Skornicki M, Parsons NB. The comparative effectiveness of a human fibroblast dermal substitute versus a dehydrated human amnion/chorion membrane allograft for the treatment of diabetic foot ulcers in a real-world setting. *Wounds.* 2017; 29(5): 125-132. 40. Frykberg RG, Gibbons GW, Walters JL, et al. A prospective, multicentre, open-label, single-arm clinical trial for treatment of chronic complex diabetic foot wounds with exposed tendon and/or bone: positive clinical outcomes of viable cryopreserved human placental membrane. *Int Wound J.* 2016; doi: 10.1111/ijw.12649. 41. Farivar BS, Toursavadkoshi S, Monahan TS, et al. Prospective study of cryopreserved placental tissue wound matrix in the management of chronic venous leg ulcers. *J Vasc Surg Venous Lymphat Disord.* 2019; 7(2): 228-233. 42. Ananian CE, Davis RD, Johnson EL, et al. Wound Closure Outcomes Suggest Clinical Equivalency Between Lyopreserved and Cryopreserved Placental Membranes Containing Viable Cells. *Adv Wound Care (New Rochelle).* 2019; <http://doi.org/10.1089/wound.2019.1028>.