

3M Science.
Applied to Life.™



3M™ Coban™ 2 and 3M™ Coban™ 2 Lite Two-Layer Compression Systems

3M has advanced the science of compression therapy by designing materials engineered with Intelligent Compression Dynamics. These materials, which comprise the Coban™ 2 and Coban™ 2 Lite Two-Layer Compression Systems, create a conformable, inelastic sleeve that supports the patient's muscle movements for effective venous return and reduction of edema.*

*Refer to instructions for use.



3M™ Coban™ 2 and 3M™ Coban™ 2 Lite Two-Layer Compression Systems

Key benefits



Easy to use: Application is fast and easy learn (For Coban™ 2 Two-Layer Compression System only)



Stays in place: Provides sustained therapeutic compression without slippage



Comfortable to wear: High patient compliance due to product's comfort during wear and low profile



Not made with natural rubber latex



Conforms to the patient's anatomy: Accommodates a variety of shapes and sizes



Provides consistent compression: Compression layer is designed to be applied at full stretch, reducing application variability



Enables a more normal lifestyle: The thin, lightweight, breathable sleeve allows patients to wear their own shoes, so they can return to their regular daily activities

Suggested applications for compression therapy

- ▶ Management of venous leg ulcers, edema, lymphedema and other clinical conditions where compression is appropriate
- ▶ Coban 2 Two-Layer Compression System (Cat #2094N) may be used for patients with an ABPI greater than or equal to 0.8
- ▶ Coban 2 Lite Two-Layer Compression System (Cat #2794N) may be used for patients who have mixed etiology, an ABPI greater than or equal to 0.5, are new to compression (where tolerance is not known), or are frail or less mobile

Ordering information

Application/Use	Name	Cat. No.	Size	Black Nylon Stocking	Rolls/box Boxes/case	Boxes/Case	HCPCS Code
Below the knee	3M™ Coban™ 2 Two-Layer Compression System	2094N	Comfort Layer: 4 in × 2.9 yd (10 cm × 2,7 m) Compression Layer: 4 in × 5.1 yd (10 cm × 4,7 m) stretched	1	2 rolls (1 of each layer) 8	6	A6441 A6452
Below the knee (ABPI ≥ 0.5)	3M™ Coban™ 2 Lite Two-Layer Compression System	2794N	Comfort Layer: 4 in × 2.9 yd (10 cm × 2,7 m) Compression Layer: 4 in × 5.1 yd (10 cm × 4,7 m) stretched	1	2 rolls (1 of each layer) 8	6	A6441 A6452
Below the knee (large circumference)	3M™ Coban™ 2 Two-Layer Compression System	2094XL	Comfort Layer: 4 in × 3.8 yd (10 cm × 3,5 m) Compression Layer: 4 in × 6.3 yd (10 cm × 5,8 m) stretched	1	2 rolls (1 of each layer) 8	4	A6441 A6452
Above the knee	3M™ Coban™ 2 Two-Layer Compression System	20096	Comfort Layer: 6 in × 3.8 yd (15 cm × 3,5 m) Compression Layer: 6 in × 4.9 yd (15 cm × 4,5 m) stretched	–	2 rolls (1 of each layer) 8	4	–
Toe boot	3M™ Coban™ 2 Two-Layer Compression System	2092	Comfort Layer: 2 in × 1.3 yd (5 cm × 1,2 m) Compression Layer: 2 in × 3 yd (5 cm × 2,7 m) stretched	–	2 rolls (1 of each layer) 8	4	–



For more information, contact your 3M Health Care Sales Representative, call the 3M Health Care Customer Helpline at 1-800-228-3957 or visit 3M.com/Coban2Layer.



3M Medical Solutions Division
2510 Conway Ave.
St. Paul, MN 55144 USA
Phone 1-800-228-3957
Web 3M.com/medical

HCPCS codes have been provided to assist you in the preparation of insurance claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.

3M and Coban are registered trademarks of 3M. Please recycle. Printed in USA. © 3M 2019. All rights reserved. 70-2011-7813-7