



Advanced
Medical
Solutions
Technology™

Skin Affix™ Topical Skin Adhesive
Designed to Optimize Wound Closure

Strong and Lasting Wound Closure

Skin Affix™ is a topical skin adhesive used to close wounds from surgical incisions, minimally invasive surgical procedures and simple, thoroughly cleansed trauma-induced lacerations.

Performance

- » Sterile 2-octyl cyanoacrylate
- » High viscosity formulation
- » Effective microbial barrier protection
- » Fast setting
- » Naturally sloughing (5-10 days)

Application

- » Versatile—foam tip allows for broad and narrow strokes
- » Promoting controlled delivery—transparent applicator
- » Easy to use—no assembly required

Patient Comfort

- » Low exothermic reaction*
- » Flexible and water-resistant


Ordering Information

Item No.	Description	Pkg.
MSC091040Z	0.4 ml	12/bx
MSC091040	0.4 ml	72/cs

Contact your Medline Representative or call 1-800-MEDLINE for more information.



Before



After one layer application



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

*Data on file.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
© 2017 Medline Industries, Inc. All rights reserved. Skin Affix is a trademark and Medline is a registered trademark of Medline Industries, Inc. Advanced Medical Solutions Technology is a trademark of Advanced Medical Solutions (Plymouth) Ltd. MKT1793632 / LIT535R / XM / XX / 3