

*The only NPWT System
powered by Hydrofiber[®] Technology*



For Optimal Incision Management

For low-to-moderate exuding wounds such as surgically closed incisions

- Powered, disposable pump with **30-day** lifetime*
- Separate pump and dressing packs
- Unique, one-way valve helps maintain pressure



Order separate pump and/or dressing packs

Dressing Size	Pad Size	Max. Incision Length	Product Code	Pack Size
12 cm x 21 cm / 4.7 in x 8.3 in	4 cm x 13 cm / 1.6 in x 5.1 in	6 cm / 2.4 in	421554	5
12 cm x 31 cm / 4.7 in x 12.2 in	4 cm x 23 cm / 1.6 in x 9.1 in	16 cm / 6.3 in	421555	5
12 cm x 41 cm / 4.7 in x 16.1 in	4 cm x 33 cm / 1.6 in x 13 in	26 cm / 10.2 in	422155	5
16 cm x 16 cm / 6.3 in x 6.3 in	8 cm x 8 cm / 3.2 in x 3.2 in		421552	5
16 cm x 21 cm / 6.3 in x 8.3 in	8 cm x 13 cm / 3.2 in x 5.1 in		421553	5
21 cm x 26 cm / 8.3 in x 10.2 in	13 cm x 18 cm / 5.1 in x 7.1 in		422156	5
26 cm x 26 cm / 10.2 in x 10.2 in	18 cm x 18 cm / 7.1 in x 7.1 in		422157	5
Pump Unit	N/A		422285	1
Pump Carrying Case	N/A		446650	1

*Battery change may be required during pump lifetime.

†As demonstrated *in-vitro*.

1. Assessment Of The In-Vitro Properties AVELLE™ Negative Pressure Wound Therapy Dressing. WHRI4520 MS128. Data on file. 2015. ConvaTec.

Call ConvaTec at **866-968-6376**

www.convatec.com