



## SC-2004-OC Sequential Circulator

Model SC-2004-OC Sequential Circulator is a gradient, sequential, pneumatic compression device, intended for the primary or adjunctive treatment of primary or secondary Lymphedema. The device is also used for the treatment of chronic venous stasis ulcers and associated venous insufficiency, as well as general treatment of swelling of the extremities. The SC-2004-OC is intended for home or hospital use.

### Key Features:

- 4-Chamber System
- Gradient (4 levels)
- Sequential
- Non-peristaltic
- Distal to proximal cycle
- Unilateral or bilateral use
- 90 seconds cycle time
- 20 - 80 mmHg pressure range
- Slim and light design
- Silent
- Easy handling
- 3 Year Warranty

### Benefits:

- Reduces Pain
- Improves Lymphatic and Vascular Flow
- Promotes Wound Healing
- Prevents Infection
- Provides Effective Prophylaxis for Venous Thrombosis
- Reduces Dynamic Edema

### Specifications:

Dimensions:  
4.5" x 11" x 7.375"  
11.43 cm x 27.94 cm x 18.73 cm  
Weight: 4.65 lb  
2.1 kg



Supply Item: E0651



### Contraindications:

Pneumatic compression is contraindicated for patients with the following:

- Congestive heart failure
- Deep Vein Thrombosis
- Inflammatory phlebitis or episodes of pulmonary embolism
- Infections in the limb, including cellulitis, without appropriate antibiotic coverage
- Active Cancer except for palliative care
- Other indications as identified by the treating physician
- Tachycardia

WARNING: Federal law restricts this device to sale by or on the order of a physician.



# Bio Compression Systems' Biocomfort Garments

Made of high quality 200 Denier Nylon Oxford coated with 3 mils of Polyurethane  
Available in a variety of sizes and styles including: over the shoulder, adjustable, custom and amputee.

Upper Extremity/4-Chamber Garment Sizing Chart								
Circumferences							Lengths	
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"
GS-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"
GS-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"
GS-3035-SH/S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"+8"
GS-3035-SH/M	Medium	27"	26"	24"	21"	18"	28"	31"+8"
GS-3035-SH/L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"+8"

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

Lower Extremity/4-Chamber Garment Sizing Chart												
Circumferences								Lengths				
Part Number	Size		Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease		
GS-3045-H	Half	Standard	N/A		29"	27"	23.5"	22"	19"	NA		
GW-3045-H		Wide			30.5"	28"	29"	23"				
GA-3045-H		Adjustable			32"	30"	23.5"	22"				
GWA-3045-H		Adj. Wide			33.5"	31"	29"	23"				
GN-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"		
GS-3045-S		Standard	31"	30"	28"	24.5"	21.5"	19.5"				
GW-3045-S		Wide	35.5"	34"	30.5"	28"	25.5"	21.5"				
GA-3045-S		Adjustable	34"	33"	31"	27.5"	21.5"	19.5"				
GWA-3045-S		Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"				
GN-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"		
GS-3045-M		Standard	32"	30"	28"	24.5"	21.5"	19.5"				
GW-3045-M		Wide	36.5"	34"	30"	28"	25.5"	22.5"				
GXW-3045		Extra Wide	44"	41"	36"	34"	30"	22.5"			29"	32"
GA-3045-M		Adjustable	35"	34"	31"	27.5"	21.5"	19.5"			28"	32"
GWA-3045-M		Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"			29"	32"
GXWA-3045	Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"			
GN-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"		
GS-3045-L		Standard	33"	32"	30"	26"	21.5"	19.5"				
GW-3045-L		Wide	37"	36"	32"	29"	24"	20.5"				
GA-3045-L		Adjustable	37"	36"	33"	29"	21.5"	19.5"				
GWA-3045-L		Adj. Wide	40"	38"	35"	32"	24"	22.5"				

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

