



SC-3008-DL Sequential Circulator Pneumatic Compression Device

The ultimate gradient sequential compression device for the treatment of Lymphedema, Venous Insufficiency, and Wound Healing. The SC-3008-DL Sequential Circulator is as unique as the patient population it services. All patients have different needs so our custom developed technology measures and then delivers pressure to each and every chamber in the garment to give the best possible treatment. The SC-3008-DL lets the user program the pressure setting in any and all of the 8 chambers, always keeping greater pressure distal so as to avoid any reflux. Complete with a compliance meter to monitor patient use, the SC-3008-DL also has a countdown clock and a timer that shuts off the pump after one hour.

Key Features:

- 8-Chamber System
- Gradient (8 levels)
- Sequential
- Non-peristaltic
- Distal to proximal cycle
- Unilateral or bilateral use
- 60 seconds cycle time
- 0 - 120 mmHg pressure range
- Slim and light design
- Silent
- Easy handling
- Timer
- Compliance meter
- Individual pressure adjustment in every chamber possible
- 3 Year Warranty

Benefits Include:

- Reduces Pain
- Improves Lymphatic and Vascular Flow
- Promotes Wound Healing
- Prevents Infection
- Provides Effective Prophylaxis for Venous Thrombosis
- Reduces Dynamic Edema

Specifications:

Dimensions: 4.5" x 11.75" x 7.75"
11.43 cm x 29.85 cm x 19.685 cm

Weight: 5.9 lb
2.68 kg

Prescription Required



Supply Item: E0652



Contraindications:

Use of this device is contraindicated for patients with the following condition:

- Congestive heart failure
- Deep Vein Thrombosis
- Inflammatory phlebitis or episodes of pulmonary embolism
- Infections in the limb, including cellulitis without appropriate antibiotic coverage
- Presence of cancer unless for palliative care



0123

Bio Compression Systems
120 West Commercial Avenue
Moonachie, NJ 07074
Phone: 800-888-0908/ Fax: 201-939-4503
email: biosystems@biocompression.com
www.biocompression.com

Bio Compression Systems' Biocomfort Garments

Made of high quality 200 Denier Nylon Oxford coated with 3 mils of Polyurethane
Available in a variety of sizes and styles including: over the shoulder,
adjustable, custom and amputee.

Upper Extremity/8-Chamber Garment Sizing Chart

		Circumferences					Lengths	
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
G8-3035-S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"
G8-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"
G8-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"
G8-3035-SH/S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"+8"
G8-3035-SH/M	Medium	27"	26"	24"	21"	18"	28"	31"+8"
G8-3035-SH/L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"+8"

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

Lower Extremity/8-Chamber Garment Sizing Chart

		Circumferences						Lengths				
Part Number	Size	Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease			
GN8-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"		
G8-3045-S		Standard	31"	30"	28"	24.5"	21.5"	19.5"				
GW8-3045-S		Wide	35.5"	34"	30.5"	28"	25.5"	21.5"				
A8-3045-S		Adjustable	34"	33"	31"	27.5"	21.5"	19.5"				
GWA8-3045-S		Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"				
GN8-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"		
G8-3045-M		Standard	32"	30"	28"	24.5"	21.5"	19.5"				
GW8-3045-M		Wide	36.5"	34"	30"	28"	25.5"	22.5"				
GXW8-3045		Extra Wide	44"	41"	36"	34"	30"	22.5"			29"	32"
A8-3045-M		Adjustable	35"	34"	31"	27.5"	21.5"	19.5"			28"	32"
GWA8-3045-M		Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"			29"	32"
GXWA8-3045	Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"			
GN8-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"		
G8-3045-L		Standard	33"	32"	30"	26"	21.5"	19.5"				
GW8-3045-L		Wide	37"	36"	32"	29"	24"	20.5"				
A8-3045-L		Adjustable	37"	36"	33"	29"	21.5"	19.5"				
GWA8-3045-L		Adj. Wide	40"	38"	35"	32"	24"	22.5"				

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb



Bio Compression Systems
120 West Commercial Avenue
Moonachie, NJ 07074
Phone: 800-888-0908/ Fax: 201-939-4503
email: biosystems@biocompression.com
www.biocompression.com